

CDpal, Version 2.16

This software is intended to be used for curve fitting of CD thermal and chemical denaturation data.

Copyright © 2017 Markus Niklasson and Patrik Lundström

This program is free software: you can redistribute it and/or modify it under the terms of the GNU General Public License as published by the free software foundation, either version 3 of the license, or (at your option) any later version.

This program is distributed in the hope that it will be useful, but WITHOUT ANY WARRANTY; without even the implied warranty of MERCHANTABILITY or FITNESS FOR A PARTICULAR PURPOSE. See the GNU General Public License for more details. You should have received a copy of the GNU General Public License along with this program. If not, see <http://www.gnu.org/licenses/>.

Programmers: Markus Niklasson and Patrik Lundström.

Date: February 21, 2017

Version: 2.16

Please report any bugs or performance issues to:

marni@ifm.liu.se Markus Niklasson

patlu@ifm.liu.se Patrik Lundström

Address correspondence to:

patlu@ifm.liu.se Patrik Lundström

Installation guide for CDpal

Windows: No installation necessary. Place folder anywhere you want and make a shortcut to the executable.

Mac: No installation necessary. Place the app in the applications directory. If you cannot start the application, Ctrl+click the application and select "Open".

Linux: Install by navigating to the CDpal directory and running the INSTALL.sh script. Afterwards, type cdpal.sh to launch the application.

Compilation information for CDpal

CDpal was compiled using Qt 5.7.0